



esoptis Sharper vision,
better results.

CLM-196 – Laboratory Colorimeter
Datasheet

CLM-196

Laboratory Colorimeter



Laboratory Colorimeter

The CLM-196 benchtop colorimeter enables spot color measurement of a wide range of production materials. This professional instrument, individually calibrated and temperature-compensated, is supplied with free application software. EOPTIS plans to continuously enhance this software by integrating new features based on specific requirements, which will be made available to all customers once development is complete.

CLM-196

Laboratory Colorimeter



Accurate, objective, daily monitoring of production.

- Practical color quality control, immediate, objective, and highly repeatable measurements;
- Sensitivity up to 60x greater than the human eye;
- Classification of measured products into batches;
- Reconstruction of the color to be produced.

CLM-196

Laboratory Colorimeter


Verify the color of raw materials and finished goods




eoptis

Measure the quality of your COLOR

Compare a color with a standard or custom database



Find the best color matching the sample



Attest your work quality with a report



Benchtop colorimeter with high reliability and long-term stability at an extremely competitive price

Conceived, developed and manufactured in ITALY



Type of materials



Cosmetics powders



Coffee



Seeds



Ceramic powders



Food liquids



Nail polishes



Fruit juices, pulps and mousses



Ice cream bases



CLM-196

Technical specifications

Color sensor	Based on the spectral value of the standard CIE 1931 color-matching functions (DIN ISO 13655 and DIN 5033)
Light source	Integrated white LEDs
Target type	Reflective (matte or glossy)
Working distance	Contact
Measurement geometry	(45°c:0°) according to CIE15:2004, ASTM E1164
Repeatability	0.03 ΔE^*_{ab} typ (30 measurements at 5s on white reference)
Illuminants	D65, D55, D50, A, C, FL2, FL7, FL11
Observers	CIE-1931 2° e CIE-1964 10°
Measurement area	Ø 10mm (standard)
Dimensions	74 mm (H) x 100mm (Ø)
Power	USB self-powered
Weight	940 g



Functional specifications

Color coordinates	CIE L*a*b, CIE XYZ, CIE L*u*v, CIE L*C*h, Yxy, sRGB, Yellowness Index, Color name.
Color differences	DeltaE*ab, dL* da* db* dC*, dX dY dZ, dL* du* dv* du*v*, dL* dC* dh, dY dx dy dxy, dR dG dB, dy
Acquisition mode	Software-triggered, single or multiple readings.
Recipe management	Ability to save parameter setups as "projects"
Report generation	Customizable reports with company logo, for analysis and history, in .pdf and .xml formats (Excel-compatible)
Database creation	Storage of measurement sequences and reference sequences
Advanced functions	Multiple measurement Average calculation of a set of measures or references Continuous measurement Creation of custom color spaces from a specific formulas

CLM-196

Accessories



Standard configuration



With «CLM-196.Tran»



With «CLM-196.Solid»

The standard configuration of the CLM-196 allows the measurement of powders, granules, and non-transparent liquids.

The CLM-196.Tran accessory consists of a dedicated eyepiece that screws onto the instrument body, together with several mechanical components that fit into the sample cup. This configuration allows the measurement of the color of transparent or semi-transparent liquids.

The CLM-196.Solid accessory consists of a dedicated eyepiece and enables the measurement of solid samples placed directly on top of the instrument.



Factory calibration and report documents



GROUP: default group												
COLOR	L*	a*	b*	SI	NO	COLOR	L*	a*	b*	SI	NO	COLOR CHECK
MSR 0020	52.85	12.55	8.44	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 24.89 dL*: -0.45 dA*: 18.36 Class: - dC*ab: -10.78 dB*: -16.80
MSR 0021	95.24	-17.20	16.00	D50	10	REF 0024	69.00	2.20	39.21	D50	10	dE*ab: 40.04 dL*: 29.25 dA*: -19.40 Class: - dC*ab: -15.77 dB*: -23.20
MSR 0017	38.58	12.13	4.54	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 31.09 dL*: -14.72 dA*: 17.95 Class: - dC*ab: -12.94 dB*: -20.69
MSR 0016	56.35	1.93	-9.00	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 35.23 dL*: 3.05 dA*: 7.74 Class: - dC*ab: -16.69 dB*: -34.24
MSR 0015	56.31	1.94	-8.98	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 35.21 dL*: 3.01 dA*: 7.76 Class: - dC*ab: -16.71 dB*: -34.21
MSR 0014	61.98	10.95	5.31	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 27.45 dL*: 6.68 dA*: 16.76 Class: - dC*ab: -13.73 dB*: -19.93
MSR 0013	62.37	12.06	17.20	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 21.60 dL*: 9.07 dA*: 17.88 Class: - dC*ab: -4.89 dB*: -8.03
MSR 0012	62.35	12.08	17.19	D50	10	REF 0027	53.30	-5.82	25.23	D50	10	dE*ab: 21.61 dL*: 9.05 dA*: 17.90 Class: - dC*ab: -4.89 dB*: -8.04
MSR 0011	77.06	-4.46	23.17	D50	10	REF 0024	69.00	2.20	39.21	D50	10	dE*ab: 19.14 dL*: 9.06 dA*: -6.66 Class: - dC*ab: -16.67 dB*: -16.84

eoptis


Factory calibration

Each device leaves our laboratories only after undergoing a complete calibration process. Thanks to our proprietary algorithms and procedures, we guarantee absolute accuracy and highly repeatable color measurements. The reliability and stability of these measurements are ensured for 12 months.

Report documents

Measurement results can be exported in PDF format, with the option to include the company logo, or in editable XML format for further analysis and processing.

«CLM19X Interface» application software




"Check Color Difference" feature

This function is used to determine the difference between two samples.

Procedure

1. Select a reference from a database or measure a reference sample
2. Set the tolerance thresholds
3. Measure the sample of interest
4. Check the difference using the graph and the tolerance level.

«CLM19X Interface» application software




«Color matching» feature

This function is used to identify which color in the database is closest to the sample of interest.

Procedure

1. Select the database to reference
2. Measure the sample of interest
3. Review the five best matches, ranked from most similar to least similar.

«CLM19X Interface» application software



«Color classification» feature

This function is used to group measurements into homogeneous color categories.

Procedure

1. Select up to 5 colors from the reference acquisitions. Each of the 5 colors defines a class
2. Measure the sample of interest.
3. Verify which class the measured sample belongs to.

Standard equipment and accessories



Standard equipment

- CLM-196 Benchtop Colorimeter
- "CLM19X Interface" application software
- Calibration certificate with RISE traceability
- White calibration target
- 1 optical glass sample cup
- Rigid, watertight carrying case with molded interior.

Accessories

- Transparent liquids kit «CLM-196.Tran»
- Solid parts kit «CLM-196.Solid»

Contact us

Eoptis Srl
Via Kufstein 15
38121 Trento (TN) Italy

+39 0461 260552

eoptis.com
info@eoptis.com



eoptis

Sharper vision,
better results.